



ABDOMINAL BINDER


ELASTIC WRAP USED TO REINFORCE THE ABDOMEN
AFTER AN ABDOMINAL SURGERY



ABD DRESSING

9 X 5 STERILE DRESSING USED TO
COVER LARGE WOUND







ACRILON PAD

FUZZY PAD TO PREVENT SKIN BREAKDOWN

ADULT DIAPERS

INCONTINENCE (URINE & STOOL)



AIR MATTRESS & MOTOR

INFLATES AND MAINTAINS AIR IN MATTRESS TO
PREVENT SKIN BREAKDOWN



ALCARE PLUS

ANTISEPTIC HAND RUB




ALCOHOL SWABS, SPONGES

SKIN PREP OR CLEANSER



**ANTI-EMBOLISM/
TED STOCKINGS**

**ELASTIC COMPRESSION
STOCKINGS USED TO
INCREASE VENOUS RETURN
& PREVENT BLOOD CLOTS**



BASIN

FOR USE IN BED BATHS










BETADINE SOLUTION SKIN CLEANSER/PREP

TOPICAL ANTISEPTIC





CASTILE SOAP

LIQUID SOAP USED FOR SOAP SUDS ENEMAS



CATHETER PLUG

USED TO PLUG UNUSED PORTS
OF 3-WAY CATH OR OTHER
MISC. TUBES






C-CLAMP

METAL CLAMP USED TO CLAMP OFF FOLEY
CATHETERS AND OTHER TUBES

CHLORA-PREP


DIFFERENT TYPE OF DISPENSER FOR
CLEANSING AGENT/ ANTIMICROBIAL AGENT





COMB

KEEPS HAIR NEAT



CONDOM CATH/TEXAS
CATH

EXTERNAL MALE CATHETER

CULTURETTE

STERILE SWAB USED TO
OBTAIN CULTURE OF WOUND



DALE CATHETER STRAP/
CATH SECURE/TAPE

CATHETER SECURING DEVICES




DENTURE CUP

PLASTIC CONTAINER TO HOLD DENTURES



DISPOSABLE CHUX

PLASTIC ABSORBANT PAD



DISPOSIBLE
PROTECTIVE EYEWEAR

UNIVERSAL PRECAUTIONS,
PROTECTS EYES FROM
FLUID SPLASHES



DISPOSABLE PULSE OX PROBE

DISPOSABLE PROBE FOR PULSE OXIMEMETRY MACHINE
USED TO OBTAIN PATIENTS OXYGEN SATURATION AND
PULSE



DISPOSABLE RAZOR

USED FOR SHAVING






DOPPLER AND GEL

USED TO AUSCULTATE BP
AND PULSES



DRAIN SPONGES

DRESSING USED WHENEVER A DRAIN IS PRESENT IN A
WOUND/INCISION



4X4 GAUZE DRESSING

STERILE DRESSING USED TO
(ALL-PURPOSE SPONGE) ABSORB OR COVER



DUODERM THICK OR THIN

PROTECTIVE DRESSING FOR SKIN BREAKDOWN






EMESIS BASIN

FOR HOLDING EMESIS


ENEMA BAG

HOLDS SOLUTION FOR
LARGE VOLUME ENEMA
FOR BOWEL CLEANSING



EPC/SCD PUMP AND CUFFS

PROVIDES INTERMITTENT PRESSURE TO
EXTREMITIES TO PREVENT BLOOD
CLOTS




EXAM GLOVES

USED WITH UNIVERSAL PRECAUTIONS



FACE SHIELD / MASK

UNIVERSAL PRECAUTIONS, PROTECTS BODY FROM
BLOOD SPLASH AND RESPIRATORY SPRAY



NDC 12462-0201-2
DIN 02239056
USP

DISPOSABLE
gent-L-tip®
PLUS
SODIUM PHOSPHATES
RECTAL SOLUTION

ready-to-use
enema

- Easy to open,
ready to use
- Prelubricated
flexible soft tip



National
Brand
Equivalent

LDT PH70166 A EXP. 02/94

Chester Labs

NET WT. 4 1/2 FL. OZ. (133ML)

REORDER NO. 101503



FLEETS ENEMA (REG.& OIL)

SMALL VOLUME ENEMA,
PREPACKAGED

Fleet.

DOCTOR
#1
RECOMMENDED

Enema
mineral oil

lubricant
laxative
ready-to-use
sodium free
latex free


Soft,
Flexible
Comfortip.



*Laxative Brand


INDWELLING FOLEY CATHETER & DRAINAGE BAG

HOLLOW TUBE TO DRAIN BLADDER, WITH BAG TO
COLLECT URINE



FORECEPS

USED TO PICK UP DRESSINGS OUT
OF WOUND, ETC.






FRACTURE PAN

SMALL BEDPAN USED FOR
PATIENT WITH FRACTURES
USUALLY WITH BEDREST RESTRICTIONS

GAIT BELT

CLOTH STRAP USED TO HELP PATIENT AMBULATE
AND MAINTAIN CONTROL



GRADUATED CYLINDER

Plastic calibrated container used
for measuring I and O



HEARING AID

ASSISTIVE DEVICE USED TO HELP HEARING



HEEL / ELBOW PROTECTORS

USED TO PROTECT BONY PROMINENCES



HEMOCULT SLIDES & SOLUTION

TEST STOOL FOR OCCULT BLOOD



HEMOSTATS

MULTI-USE INSTRUMENT FOR CLAMPING AND HOLDING



HEMOVAC

CLOSED WOUND SUCTION




ICE PACK

REFILLABLE BAG FOR ICE COMPRESS




INCENTIVE SPIROMETER

DEVICE TO EXPAND LUNGS WITH
INHALATION, PREVENTS LAZY
LUNGS, ATALACTASIS



IRRIGATION SET

PLASTIC CONTAINER WITH SYRINGE
FOR WOUND IRRIGATION



Directions for Use:

1. Remove seal cap.
2. Remove syringe from solution container. Fill container with clean irrigation solution and insert syringe into solution container.
3. Place syringe tip on wound, plastic side down.
4. Place solution tray on solution side.
5. Gradually lift tray from the patient.
6. Discontinue irrigation when solution is exhausted.

Contents:

1. Disinfect solution container
1. Disinfect plastic injection bulb
1. Irrigation container
1. Piston syringe
1. Low Adherent Puff Filter Culture
1. Syringe
1. Post and (optional) tubing
1. Solution

Caution:

Do not reuse.

Caution: Contents of tray are sterile unless package has been opened or damaged.

Caution: Federal (21 CFR) law restricts this device to sale by or on the order of a physician.

© Bard, Inc.
Charleston, SC 29404
1-800-328-4444
www.bard.com

LOT: NGSK0970
EXP. BY: 2013-06

Part No. 750307
Formerly 0508480






JACKSON PRATT DRAIN

SOFT TUBE THAT IS PLACED IN AN OPERATIVE SITE TO DRAIN BLOOD AND FLUID FOLLOWING SURGERY. THE TUBE IS ATTACHED TO A SMALL COMPRESSED BULB FOR SUCTION.



KERLEX FLUFF SPONGES


USED FOR WET/DRY DRESSING, ABSORPTION
AND PACKING ALSO FOR
PACKING



KERLIX ROLL

STERILE GAUZE FOR WRAPPING
AND PACKING WOUNDS





K-PAD OR AQUA-K

CIRCULATES WARM WATER, DRY HEAT